Characterization of Clearance Mechanisms in an

All-Human Cell Based Tri-Culture System

40-

30-



Stephanie Piekos¹, Jessica R. Weaver², Cody Thomas¹, Alexander Byer-Alcorace¹, Justin J. Odanga², Kristina K. Wolf³, Jingsong Chen², Jung Bok Lee², Edward L. LeCluyse³, and Mitchell E. Taub¹ ¹Non-clinical DMPK, Boehringer Ingelheim Pharmaceuticals Inc., Ridgefield, CT

²Institute of Regenerative Medicine, LifeNet Health, Virginia Beach, VA


³Research and Development, LifeSciences Division, LifeNet Health, Research Triangle Park, NC





Abstract


Routine use of primary human hepatocytes (PHH) for clearance estimations of low-turnover compounds has been restricted in part by the lack of a suitable technology that is convenient, supports a multitude of donor lots, and maintains the metabolic pathways of PHHs over prolonged culture periods. To address these limitations, an all-human tri-culture system has been developed comprised of cryopreserved PHHs and primary feeder cells (FCs). FCs are stromal and endothelial cells at a 1:1 ratio. Frozen FCs were thawed and seeded onto 24-well collagen coated plates and allowed to attach. Cryopreserved adult PHHs were then thawed and plated onto the FCs, creating a Tri-Culture System (TCS). Basic morphology and functionality, including formation of functional bile canaliculi, albumin (Alb) synthesis, and urea production were evaluated over a 4-week culture period. Cytochrome P450 (CYP), UDPglucuronosyltransferase (UGT) and aldehyde oxidase (AO) activities were also determined over a 2-4 week period. The ability of the TCS to predict in vivo non-renal clearance of two slowly cleared compounds was also evaluated. PHHs from multiple adult donor lots in the TCS exhibited a healthy and stable morphology, including multicellular cluster formation, for up to 30 days in the 24-well format. Extensive anastomosing networks of bile canaliculi were identified with CDFDA staining after 5 days in culture with well-formed tight and gap junctions that remained stable throughout the rest of the culture period. Albumin and urea production levels were maintained in the TCS over a 4-week period. CYP1A2, CYP2C9, CYP2D6, and CYP3A4 specific activities were maintained at stable levels over the culture period. E2-3G production (UGT1A1) was stable between day 5 and 28 of culture, while AZT-G formation (UGT2B7) was well-maintained but exhibited greater variability over the 4-week culture period. AO specific activity was greater in the TCS than in suspended PHHs and was stable over at least a 14-day period. The TCS was also able to predict the hepatic clearance of alprazolam and tolbutamide within 2-fold of reported in vivo values. Overall, these results show that the TCS represents a convenient, stable, all-human hepatic culture system that maintains both hepatocellular morphology and various metabolic functions for up to 1 month.


Results Hepatocytes in the TCS retain functionality for up to four weeks in culture Albumin Urea **day** 30-












96 120 144 168 192

72

Time (h)



- 4. Assess inducibility of major CYP enzymes in response to prototypical inducers by measuring gene expression
- 5. Evaluate the ability of the TCS to predict in vivo clearance of slowly metabolized compounds by incubating with tolbutamide and alprazolam for seven days

¹Chan TS et al. Drug Metab Dispos. 2013 Dec; 41(12).

Conclusions

120 144 168 192

96

Time (h)

72

- The TCS maintains hepatic function and high levels of drug metabolizing enzyme activity for up to four weeks in culture
- The activity levels of AO, UGT1A1, and UGT2B7 are significantly higher in TCS than in suspended hepatocytes of the same donor
- Expression levels of CYP enzymes in the TCS are highly inducible in response to prototypical inducers
- The TCS may be a useful in vitro model to predict in vivo clearance of slowly metabolized compounds