

# Accurately Predict *In Vivo* Hepatic Clearance with TruVIVO<sup>®</sup>

Accurate *in vitro* prediction of hepatic clearance ( $Cl_{int}$ ) is pivotal for successful drug development. Designing drugs for oral administration with low first-pass metabolism and longer half-lives presents challenges in predicting slow metabolism and low  $Cl_{int}$  values. TruVivo, a pioneering all-human 2D+ hepatic system, provides a solution for obtaining more relevant  $Cl_{int}$  data with high IVIVE accuracy. Features of TruVivo include:

- **All-Human Cells:** TruVivo combines pre-qualified primary human hepatocytes with human feeder cells of stromal and endothelial origin
- **Human Relevance:** Hepatocytes in this system maintain their morphology, metabolic health, drug metabolizing capacity (phase I and phase II), and transporter function for more than 14 days
- **Broad Applicability:** TruVivo can provide clearance data for high, intermediate, and low  $Cl_{int}$  compounds in a single model

## Benefits of Our Services:



**Accurate** and **reliable** data for informed decision-making



**Fast turnaround times** to keep your drug development on track




**Unsurpassed expertise** and quality control from the donor gift to assay testing



**Collaborative** approach to ensure the study addresses your research questions

## Additional LifeNet Health LifeSciences Services | [Click to learn more:](#)


- *In vitro* assay services
- DMPK/ADMET assays
- OECD safety studies
- Cytotoxicity assays



## Pre-Study Evaluation Study Design

ASSAY PARAMETER	PROTOCOL
<b>Test system</b> (panel of donors for client selection)	TruVivo, all-human 2D+ hepatic system
<b>Qualified hepatocyte lots</b>	Metabolism, genotyping, clinical background
<b>Plate format</b>	24-well
<b>Replicates</b>	3
<b>Cell density</b>	0.3 x 10 <sup>6</sup> /well
<b>Test Cpd concentration</b> (can be adjusted)	1 $\mu$ M
<b>Standard sample times</b> (low, medium, and high clearance)	0, 1, 2, 4, 8, 22, 30, 48, 72, 96hr (up to 7 days is possible)
<b>Amount of test compound required</b>	2mg
<b>Cell health markers</b>	Urea production and LDH leakage
<b>Bioanalytical</b>	LC/MS/MS
<b>Reference Cpd<sub>s</sub> (Cl<sub>int</sub>)</b>	2 (high and low)
<b>Time to complete</b>	3-4 weeks
<b>Regulatory</b>	Non-GLP or GLP compliant
<b>Deliverables</b>	Full report, Cl <sub>int</sub> value, T <sub>1/2</sub> , with graphs and tables where appropriate

## Sample Plate Layout



Example 24-well plate layout for six test compounds and two reference compounds

## Potential Reference Compounds

Compound name	Cl <sub>int</sub> category
Diazepam	Low
Disopyramide	Low
S-Warfarin	Low
Metoprolol	Low
Naloxone	High
Bupirone	High
Verapamil	High

## How can we help?

Talk with one of our experts for help with general inquiries, protocol details, or becoming a new client.

[LNHLifesciences.org/information](https://lnhlifesciences.org/information)  
1-888-847-7831 (US & Canada)  
[cells\\_tissues@lifenethealth.org](mailto:cells_tissues@lifenethealth.org)