


In Vitro Tri-culture Hepatic Model to Evaluate Human Relevance of Chemical-Induced Thyroid Toxicity

Ahtasham Raza¹, Kristina K. Wolf², Mercedes Biven³, Tammy Stone², Stephanie Kellum¹, Edward LeCluyse², Jessica LaRocca³, Raja Settivari¹ and Shadia Catalano¹.
 Corvea Agriscience, Newark (DE)¹. Corvea Agriscience, Indianapolis (IN)³. LifeNet Health, Durham (NC)².

INTRODUCTION

Background

- ECHA and EFSA published a guidance for TH disruptors in animals should be investigated for human relevance to enhance the safety measure of novel plant protection products.
- Advancements in NAMs provides *in vitro* test systems to evaluate the catabolism rate of TH. It has the potential to prediction species differences response, gender disparity and population variability.




Adverse outcome pathway (AOP) of Nuclear receptor agonist inducing thyroid perturbation. ME (green) key events (blue) and associated events (grey) are the network of TH signaling that involves feedback interactions of the hypothalamic-pituitary-thyroid (HPT) axis, circulatory system, liver, and other tissues Noyes et al. (2019).

AIMS

- The model was then used to evaluate TH perturbation (key and associated events) after assault of known nuclear receptor agonists.


METHODS

- Repeated 7-day exposure to Phenobarbital (PB, 100 and 500µM), Rifampicin (RIF,10µM), Pregnenolone 16α-carbonitrile (PCN, 20µM), Polychlorinated biphenyl 153 (PCB153, 30µM).
- On Day 9 of culture, addition of labeled T4 (0.05 µM, and 0.1 µM rat and human respectively)
- 3 Human (2 Male,1 Female) and 1 Rat (SD, male).
- LC-MS/MS analysis of T4, T4G, T4S, T3 (in media).
- T4 Intrinsic clearance (CL_{int} (*in vitro*) calculation (Richardson et al. (2014).
- qRT-PCR for phase I and II enzymes.




RESULTS

T4 METABOLISM RAT vs HUMAN




T4 metabolism in Rat and Human: T4G production in both species showed a constant increased of T4G over 24-hours. Significant higher rate of production was seen in rat (SD) compared to human. PCB153 and PCN exposure has further increased the production in rat over 24-hours. While in human only PCB153 exposure has increased T4G production (relative to DMSO control). T4S no peak was observed in both species. T3 production in both species was at a similar rate with a consistent increase over 24-hours. RIF exposure has increased the production in human after 24-hours (not seen in rat). T4 concentration in media decreases after exposure to PCB153 (11%), PCN (44%), and PB500 (75%) after 24-hours in rat. While in human, the decrease was noticed after exposure to PB500 (63.5%), RIF (75%) and PCB153 (64%). T4 clearance rate Intrinsic clearance (CL_{int} (*in vitro*) of T4 in the rat model after PCB153 and PCN treatment was calculated at 4.2 and 1.9 µL/min/10⁶ cells, respectively. While in human tri-cultures, Intrinsic clearance (CL_{int} (*in vitro*) of T4 after PB and RIF was calculated at 1.2 and 1.4 µL/min/10⁶ cells, respectively.

Cyp mRNA EXPRESSION RAT vs HUMAN



Ugt mRNA EXPRESSION RAT vs HUMAN



CONCLUSION

- PCB153 MoA:** Increased T4G production in both rat and human over 24-hours, T3 production was unaffected in both species. Higher T4 clearance rate in rat than in human over 24-hours. Cyp2B & 3A upregulated in both species. Ugt2B upregulated in rat.
- PCN MoA:** Increased T4G production in rat but not in human over 24-hours. T3 production slightly decreased in rat (not in human). Higher T4 clearance rate in rat (not in human). Cyp3A and Ugt2B1 expression upregulated in rat only (no Cyp/Ugt upregulation in human).
- RIF MoA:** No effect on T4G production in both species, but a significant T3 increase was noticed in human. Slightly higher T4 clearance rate was noticed in human only. Cyp2B/3A and Ugt1A1 genes upregulated in human only.

Compounds	Rat	Human
	CL _{int} in vitro (µL/min/10 ⁶ cells)	CL _{int} in vitro (µL/min/10 ⁶ cells)
Vehicle	0.630	0.840
PB-100	0.840	1.050
PB-500	0.840	1.261
CITCO	0.630	1.050
DEX	1.471	1.050
RIF	0.840	1.471
PCN	1.261	1.050
PCB153	4.202	1.050

Intrinsic clearance (CL_{int} (*in vitro*) of T4 in the rat and human Tri-culture model

Species	Compound	Nuclear receptor agonist	T4 Clearance											
			T4G	T3	T4	CYP1A1	CYP2B1	CYP3A23	UGT1A1	UGT1A6	UGT2B1			
Rat	PB	CAR (human & rat)	-	-	-	-	-	-	-	-	-	-	-	-
	CITCO	CAR (human)	-	-	-	-	-	-	-	-	-	-	-	-
	DEX	PXR (human & rat)	-	D	-	-	-	-	-	-	-	-	-	-
	RIF	PXR (human)	-	D	-	-	-	-	-	-	-	-	-	-
	PCN	PXR (rat)	-	D	-	-	-	-	-	-	-	-	-	-
Human	PCB153	CAR/PXR (human & rat)	I	D	I	-	-	-	-	-	-	-	-	-
	Compound	Nuclear receptor agonist	T4G	T3	T4	CYP1A1	CYP2B1	CYP3A23	UGT1A1	UGT1A6	UGT2B1			
	PB	CAR (human & rat)	-	-	-	-	-	-	-	-	-			
CITCO	CAR (human)	-	-	-	-	-	-	-	-	-				
DEX	PXR (human & rat)	-	D	-	-	-	-	-	-	-				
RIF	PXR (human)	-	D	-	-	-	-	-	-	-				
PCN	PXR (rat)	-	D	-	-	-	-	-	-	-				
PCB153	CAR/PXR (human & rat)	I	D	I	-	-	-	-	-	-				

SUMMARY

- T4G production rate was increased over 24-hours in both species. The production rate was much higher in rat than in human.
- T3S was not observed and requires further investigation.
- Species difference in T4 clearance was observed with PCN (rat-sensitive), PCB153 (rat-sensitive) and RIF (human-sensitive).
- PB, CITCO, DEX did not show any response in T4G. Though no changes in T4-clearance in either species, warranting further investigation and technical troubleshooting.
- Consistent Cyp and Ugt gene expression data of known NR agonist among both human and rat species.

Noyes et al. (2019). Evaluating chemicals for thyroid disruption: Opportunities and challenges with *In Vitro* testing and adverse outcome pathway approaches. *Environ Health Perspect*. 127(9). Richardson et al. (2014). *In vitro* metabolism of thyroxine by rat and human hepatocytes. *Xenobiotica*, 44(5), 391-403.